Germany  
**Feller**


BAP 250 D

**Electronic  
Personal  
Scale**


## BAP 250 D

Before operating this unit, Please read the instruction manual carefully.  
This product is for household use only.


## DESCRIPTION

**Note:**

1. Please take out the read-out display when setting time and alarm
2. Please do not cover the infrared emitter and receiver when weighing under wireless transmission mode.

1. Read-out display
2. Infrared emitter
3. Main body


1. Real time/setting
2. Alarm/increase
3. Infrared receiver

4. To weigh/decrease
5. Recall/Compare/Start test

**Note:**

1. Please take out the read-out display when setting time and alarm.
2. Please do not cover the infrared emitter and receiver when weighing under wireless transmission mode.

- |                     |                              |
|---------------------|------------------------------|
| 1. Read-out display | 5. Real time/setting         |
| 2. Infrared emitter | 6. To weigh/decrease         |
| 3. Main body        | 7. Recall/Compare/Start test |
| 4. Alarm/increase   |                              |

**Note:**

1. Please take out the read-out display when setting time and alarm
2. Please do not cover the infrared emitter and receiver when weighing under wireless transmission mode.

**Dear Customer,**

Welcome to use Feller Body Fat/Hydration Monitor Scale. This product is specially designed to test body fat/hydration/muscle/bone, help keep a close eye on the daily change of your body fat, hydration, muscle, bone and weight, providing information for your reference. Base on the analysis of your height, weight and age etc., it tells your daily energy requirement, and helps control your weight. Please read this instruction carefully before use.

**USER-FRIENDLY TIPS**

1. During Measuring
2. The device is designed to allow auto step on. Please always ignore the first reading and only start taking effective record from the second time weighing.
3. Always remove your shoes and socks, wipe clean your feet before using.
4. Step on the weighing platform gently.
5. It is advisable to take measurement at the same time of the day.
6. Measured result could be misleading after intensive exercise, excessive dieting or under extreme dehydration condition.
7. Always weigh / use the scale on a hard and flat surface. When you are measuring, please keep still.
8. The measured data of the following people may have deviation:
  - Children under 7 or adult over 99 (can only use the device under normal weighing mode)
  - Adult over 70
  - Body building athletes or other occupational athletes.
9. People with the following conditions should not use the device:
  - Pregnant women
  - Those with symptom of edema
  - Those on dialysis treatment
  - Those who use heart pacemaker or with other implanted medical device.

**SAFETY WARNING**

1. The Fat% Hydration% content evaluated by the Fat Monitor Scale is only for reference (not for medical purpose). If your Fat% Hydration% content is over/below normal level, please consult your doctor for more advice.
2. Keep away from strong electro-magnetic field when using the scale.
3. It could be slippery if you weigh on a wet surface. Moreover, if the weighing platform surface is glazed, never step on the edge of the scale and do keep balance during your weighing.

**FEATURES / SPECIFICATIONS**



1. Infrared transmission
2. Dual display feature: wired red-out display or wireless read-out display

3. Auto return to real time clock if no operation
4. Multifunction: test body fat %, body hydration %, body muscle %, body bone% and suggest calorie intake requirement.
5. Five activity levels selectable
6. Built-in memorizer for 12 family users and 1 guest user
7. Optional normal weighing mode (bypassing parameter setting): Auto step on weighing and consecutive add-on weighing
8. Low battery and overload indication.
9. Last immediate memory recall; compare current weight reading to the last record
10. Keypad tone

### PERFORMANCE INDEX


Item Index	Fat%	Hydration %	Muscle %	Bone %	Calorie
Index Range	4.0% ~ 60.0%	27.5% ~ 66.0%	20.0% ~ 56.0%	2.0% ~ 20.0%	
Division	0.1%	0.1%	0.1%	0.1%	1KCAL

### LCD DISPLAY


### TO INSTALL BATTERY

1. Powered by 4 x 1.5V AAA batteries. 2 for scale, 2 for read-out display
2. Take out used batteries, with the aid of a sharp object if needed. Install new batteries by putting one side of the battery down first and press down the other side. (Please take note to the polarity)


### TO SET REAL TIME CLOCK

When scale displays real time, press "SET" for 3 seconds to set time.


### TO SET ALARM

When scale displays real time, press “🔔” to start setting alarm.

### TO SET WEIGHT UNIT (UNDER NORMAL WEIGHING MODE)

If your scale is available (or applicable) in kg/lb/st, you can set to the one you prefer as follow.


After selecting the preferred unit, the scale will be activated when the scale is turned on again.

- When setting lb/st as weight unit, the default height unit is feet/inch;
- When setting kg as weight unit, the default height unit is cm.

Remark: The unit system of scale depends on sales destination or the requirement of distributor.


### TO WEIGH

#### 1. Wired transmission mode


#### 2. Wireless transmission mode


A. To weigh under wireless transmission mode.

**B.** You can weigh with read-out display adsorbed to metal surface.


**C.** You can weigh with read-out display hanged on a nail on the wall.

### CONSECUTIVE ADD-ON WEIGHING

After digits stabilize, more additional loadings can be conducted. Ensure each add-on loading has to exceed 2kg to get a new reading. Example: If you want to get a weighing value of your baby, you can weigh yourself first, and then hold your baby together, without stepping off the scale. Then you will get the total weight reading. Subtracting the previous reading of yourself from the new reading, you will get the difference as your baby's weighing value.

### FAT / HYDRATION / MUSCLE / BONE / CALORIE MEASURING MODE

#### (1) Setting personal parameters

1. At the first time when you install battery and switch on, all default parameter are set as follow:


Gender	Height	Age	Activity level
Female	165cm	25	

2. Parameters can be set in the following range:

Memory	Gender	Height	Age	Activity level
1 ~ 12 and VIP	male/female	100 - 250cm 3'03.5" ~ 8'02.5"	7 ~ 99 Year	

3. Parameter setting

**Example:** the 5th user, female, 166cm, 26 years old can set her parameters as follow:

### Note:

- When setting parameters, please take out the read-out display
- During the parameter-setting, you can always tap on the weighing platform any time to set the scale to zero-point. You can start measuring based on the parameters just set/updated into the scale.
- It has 12 user codes for family users, and 1 user code for guest user. The personal data of guest user cannot be memorized.
- The setting range of age should be from 7 to 99. For child aged from 7 to 17, LCD will display the icon of Child. For adult aged from 18 to 99, LCD will display the icon of Adult.

### • Activity level:

Sedentary/very inactive: little or no exercise

Moderate activity: always stand or walk

Light activity: Mental workers

Very active: manual workers


Extremely active: athletes


## 2. To Start Measuring

To ensure accuracy, please remove your shoes & socks before stepping on. Align two feet respectively onto the two sides of scale electrodes. DO NO step off before LCD displays test result.

**Example:** the 5th user, female, 166cm, 26 years old can measure as follow:

TO RECALL TEST RECORD


**Remark:** no record for guest user

### TO COMPARE TEST RESULT (FOR FAMILY USERS ONLY)

After body fat measuring, LCD shows current test result. Meanwhile you can press START to view current test result and last test result. Last result will be displayed in the real time display zone and flash. Press ▲ or ▼ to view current and last body fat, body hydration, body bone, and body muscle result respectively. Press SET to return real time display.

### WARNING INDICATIONS

#### Low Battery power for read out display

Please replace battery.

#### Low battery power for main unit.

Please replace battery.

#### Over- load indication

The weighting subject on the platform exceeds the maximum capability of scale.

Please stepp off to avoid damage.

#### High Fat% Indication

The Fat% is too high. Please watch your diet and do more exercise.

#### Re-test

Error exists, please re-test to get the correct result.

## ABOUT FAT%

### 1. About Fat%

- Fat% is an index % of fat content in human body.
- The fat content in human body is of a certain level, overabundance of fat will bring bad effects on metabolism, thereby, various illnesses may invade. Here goes an old saying, "Obesity is hotbed of illness".

### 2. Working Principles

As electrical signal of certain frequency transmits through human body, the impedance of fat is stronger than that of muscle or other human tissues. Through the use of biological resistance analysis, this product will pass a safe frequency of low voltage electrical current through human body, thus, the Fat% will be reflected according to the resistance differences when the electrical signal transmitting through human body.

### 3. Fat / Hydration Content Fitness Assessment Chart

After your weighing, you can use the following Assessment Chart for reference:

Status	Too lean	Lean	Normal	Fat	Too Fat
Key Line Indication	.....i	.....iiii	.....iiiiiiii	.....iiiiiiiiiiii	.....iiiiiiiiiiiiiiii

## ABOUT CALORIE

### 1. Calorie & Energy Consumption

Calorie: Energy Unit

The data reading indicate the energy you need to sustain the daily metabolism and ordinary activities required for your weight, height, age and sex.

### 2. Calorie & Losing Weight

Control the calorie intake combined with proper exercise is an effective way to lose weight. This is a very simple theory. As the daily energy intake cannot support the energy consumption of the human body, the saccharide and fat stored inside will decompose and provide energy to human body. That leads to weight loss.

## CLEANING & MAINTENANCE

1. Use a slightly damp soft fabric to wipe clean the scale. DO NOT immerse it into water. Any use of chemical or corrosive detergent is not acceptable.
2. If the plastic components of the scale contacts with grease, spicy sauce, vinegar or other irritating substance, please clean it at once (Please keep the scale away from any fruit juice).
3. Always weight / use the scale on a hard and flat surface. DO NOT use on carpets or mats.

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